

BPC BLUE Premium (1L) DUO Components List

ART NO 24-0500-02



The BPC Blue Premium (1L) DUO is composed of:

Incubation unit:

20 GL45 1L glass bottles (ART NO 95-0001-00)

Thermostatic water bath III (ART NO 21-0104-01)

Water bath lid and base tray (ART NO 22-0105-01, 21-0107-01)

Motor controller unit and its power adaptor, signal cable, and power splitter. (ART NO 01-0109-0, 17-0002-01, 94-0021-00, 01-0311-01)

19 Brushless DC Motor cables (15x250 mm, 3x500 mm, 1x1500mm) (ART NO 94-0018-00, 01-0103-04, 01-0104-02)

18 Brushless DC Motors (ART NO 01-0102-02)

20 Axis couplings (ART NO 01-0112-01)

20 Stirrers for GL45 (1L) (ART NO 21-0103-03)

20 Stirrers for GL45 (GAA) (1L) (ART NO 24-0102-02)

CO2 Absorption Unit:

2 Bottle holders (9x250 mL) (ART NO 01-0204-03)

20 Glass bottles (250 mL) (ART NO 12-0101-01)

20 Bottle nuts GL45 (ART NO 01-0202-02)

20 Lids GL45 (ART NO 01-0203-01)

20 Gas absorption attachments (GAA) (ART NO 24-0201-01)

20 GAA Funnel (MK2) (ART NO 24-0202-01)

20 GAA Condensation traps

BPC Core Unit

2 BPC Core, power adapters and 3 ethernet cables. (ART NO 28-0301-01)

40 Flow cell units 2mL (ART NO 12-0302-01)

40 Flow cell units 9mL (ART NO 01-0305-02)

FCU volume sheet for calibration

1 Ethernet desktop switch (ART NO 11-0008-01)

Additional Components

2 Festo tubing 50m (ART NO 08-0610-01)

20 Tubing stoppers (ART NO 21-0406-02)

20 Push-in valves (ART NO 24-0401-01)

160 Check valves (ART NO 90-1009-01)

80 Push-in connectors 6mm (ART NO 24-0402-02)

3 Gas bags (ART NO 105-0003-00)

1 Gas distribution manifold (3) (ART NO 103-0001-00)

3 Gas distribution manifold (6) (ART NO 12-0202-06)

12 Soft binders 7/180 mm (ART NO 21-0106-02)

40 Multi-coloured marker clamps (ART NO 21-0405-02)

Funnel (ART NO 21-0407-01)

Number marking stickers & Syringe

Bottle/tube opening tool (ART NO 24-0403-01)