

BPC Blue Anaerobic Light (2L) DUO Components List

ART NO 25-0500-01



The BPC Blue Anaerobic Light (2L) DUO is composed of:

Incubation unit:

- 10 GL45 2L glass bottles (ART NO 102-0003-00)
- Thermostatic water bath III (ART NO 21-0104-01)
- Water bath lid and base tray (ART NO 21-0105-02, 21-0107-01)
- Motor controller unit and its power adaptor, signal cable, and power splitter. (ART NO 01-0109-0, 17-0002-01, 94-0021-00, 01-0311-01)
- 10 Brushless DC Motor cables (6x250 mm, 3x500 mm, 1x1500mm) (ART NO 94-0018-00, 01-0103-04, 01-0104-02)
- 9 Brushless DC Motors (ART NO 01-0102-02)
- 10 Axis couplings (ART NO 01-0112-01)
- 10 Stirrers for GL45 (2L) (ART NO 02-0102-02)
- 10 Lids GL45 (ART NO 01-0203-01)

BPC Core Unit

- 1 BPC Core, power adaptor and 3 ethernet cables. (ART NO 28-0301-01)
- 20 Flow cell units, 2 mL or 9mL (ART NO 12-0302-01 or 01-0305-02)
- FCU volume sheet for calibration
- 1 Ethernet desktop switch (ART NO 11-0008-01)

Additional Components

- 1 Syringe
- 1 Festo tubing 50m (ART NO 08-0610-01)
- 20 Tubing stoppers (ART NO 21-0406-02)
- 10 Push-in valves (ART NO 24-0401-01)
- 40 Check valves (ART NO 90-1009-01)
- 40 Push-in connectors 6mm (ART NO 24-0402-02)
- 6 Soft binders 7/180 mm (ART NO 21-0106-02)
- 20 Multi-coloured marker clamps (ART NO 21-0405-02)
- Funnel (ART NO 21-0407-01)
- Number marking stickers
- Bottle/tube opening tool (ART NO 24-0403-01)