

# BPC Blue Premium Light (2L) Components List

ART NO 25-0000-03



The BPC Blue Premium Light (2L) is composed of:

## Incubation unit:

- 10 GL45 2L glass bottles (ART NO 102-0003-00)
- Thermostatic water bath III (ART NO 21-0104-01)
- Water bath lid and base tray (ART NO 21-0105-02, 21-0107-01)
- Motor controller unit and its power adaptor, signal cable, and power splitter. (ART NO 01-0109-0, 17-0002-01, 94-0021-00, 01-0311-01)
- 10 Brushless DC Motor cables (6x250 mm, 3x500 mm, 1x1500mm) (ART NO 94-0018-00, 01-0103-04, 01-0104-02)
- 9 Brushless DC Motors (ART NO 01-0102-02)
- 10 Axis couplings (ART NO 01-0112-01)
- 10 Stirrers for GL45 (2L) (ART NO 02-0102-02)
- 10 Stirrers for GL45 (GAA) (2L) (ART NO 24-0103-02)
- 10 Bottle nuts GL45 (ART NO 01-0202-02)
- 10 Lids GL45 (ART NO 01-0203-01)

## CO<sub>2</sub> Absorption Unit:

- 1 Bottle holders (9x250 mL) (ART NO 01-0204-03)
- 10 Glass bottles (250 mL) (ART NO 12-0101-01)
- 10 Gas absorption attachments (GAA) (ART NO 24-0201-01)
- 10 GAA Funnel (MK2) (ART NO 24-0202-01)
- 10 GAA Condensation traps (ART NO 24-0203-02)

## BPC Core Unit

- 1 BPC Core Light, power adaptor and ethernet cable. (ART NO 29-0301-01)
- 10 Flow cell units 2mL (ART NO 12-0302-01)
- 10 Flow cell units 9mL (ART NO 01-0305-02)
- FCU volume sheet for calibration

## Additional Components

- 1 Syringe and Number marking stickers
- 1 Festo tubing 50m (ART NO 08-0610-01)
- 10 Tubing stoppers (ART NO 21-0406-02)
- 10 Push-in valves (ART NO 24-0401-01)
- 40 Check valves (ART NO 90-1009-01)
- 20 Push-in connectors 6mm (ART NO 24-0402-02)
- 3 Gas bags (ART NO 105-0003-00)
- 4 Gas distribution manifold (3) (ART NO 103-0001-00)
- 6 Soft binders 7/180 mm (ART NO 21-0106-02)
- 20 Multi-coloured marker clamps (ART NO 21-0405-02)
- Funnel (ART NO 21-0407-01)
- Bottle/tube opening tool (ART NO 24-0403-01)