

# GE III Standard (1L) DUO Components List

ART NO 28-0500-01



The Gas Endeavour III Standard (1L) DUO is composed of:

## Incubation unit:

20 GL45 1L glass bottles (ART NO 95-0001-00)

Thermostatic water bath III (ART NO 21-0104-01)

Water bath lid and base tray (ART NO 22-0105-01, 21-0107-01)

Motor controller unit and its power adaptor, signal cable, and power splitter. (ART NO 01-0109-0, 17-0002-01, 94-0021-00, 01-0311-01)

19 Brushless DC Motor cables (15x250 mm, 3x500 mm, 1x1500mm)  
(ART NO 94-0018-00, 01-0103-04, 01-0104-02)

18 Brushless DC Motors (ART NO 01-0102-02)

20 Axis couplings (ART NO 01-0112-01)

20 Stirrers for GL45 (1L) (ART NO 21-0103-03)

## CO2 Absorption Unit:

2 Bottle holders (9x250 mL) (ART NO 01-0204-03)

20 Glass bottles (250 mL) (ART NO 12-0101-01)

20 Bottle nuts GL45 (ART NO 01-0202-02)

20 Lids GL45 (ART NO 01-0203-01)

## BPC Core Unit

2 BPC Core, power adapters and 3 ethernet cables. (ART NO 28-0301-01)

40 Flow cell units 2mL (ART NO 12-0302-01)

40 Flow cell units 9mL (ART NO 01-0305-02)

FCU volume sheet for calibration

1 Ethernet desktop switch (ART NO 11-0008-01)

## Additional Components

1 Syringe

1 Festo tubing 50m (ART NO 08-0610-01)

20 Tubing stoppers (ART NO 21-0406-02)

20 Push-in valves (ART NO 24-0401-01)

160 Check valves (ART NO 90-1009-01)

80 Push-in connectors 6mm (ART NO 24-0402-02)

12 Soft binders 7/180 mm (ART NO 21-0106-02)

40 Multi-coloured marker clamps (ART NO 21-0405-02)

Funnel (ART NO 21-0407-01)

Number marking stickers

Bottle/tube opening tool (ART NO 24-0403-01)